

Recombinant Platelet/Endothelial Cell Adhesion Molecule (PECAM1)

Catalog No.: TP08941

50µg

Sequence Information

Species: Rat

Gene ID:29583

Swiss Prot:Q3SWT0

Synonyms:CD31; EndoCAM; GPIIA'; PECA1

Residues:Arg459~Glu624

DPVTFDTDKPTRDMEYQCIVDNCHSHPEVRSEILRVKVIAPVDEVTISILSGNDV
QSGDEMVLRCVKEGTGPVTFQFYKEKEGRPFHEETVNDTQVFWHHEQTSKEQE
GQYYCTAFNRASIVTSLRSGPLTVRVFLAPWKKGLIAVWVIGVIAALIVAAY
YFLRKAKAKQKPVE

Product Information

Source: Prokaryotic expression.

Host: *E. coli*

Tags:N-terminal His and GST Tag

Subcellular Location: Membrane; Single-pass type I membrane protein. Cell junction.

Purity: >95%

Traits: Freeze-dried powder

Buffer formulation:PBS, pH7.4, containing 0.1% SKL, 5% Trehalose.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 6.7

Predicted Molecular Mass: 51.2kDa

Accurate Molecular Mass: 51kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH₂O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is

less than 5% within the expiration date under appropriate storage condition.

[IDENTIFICATION]




Figure 1. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for in vitro and research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.