

## Recombinant Secondary Lymphoid Tissue Chemokine (SLC)

Catalog No.: TP01459

50µg

### Sequence Information

**Species:** Human

**Gene ID:**6366

**Swiss Prot:**O00585

**Synonyms:**C-C motif chemokine 21, 6Ckine,  
CKb9, ECL, SCYA21, SLC, TCA4.

**Residues:**Ser24~Pro134

SDGGAQDCCLKYSQRKIPAKVVRSYRKQEPSLGCSIPAILFLPRKRSQAELCAD

PKELWVQQLMQHLDKTPSPQKPAQGCRKDRGASKTGKKGKGSKGCKRTERSQTP

KGP

### Product Information

**Source:** Recombinant expression.

**Host:** *E.coli*

**Tags:** N-terminal His-Tag

**Subcellular Location:** Secreted.

**Purity:** >90%

**Traits:** Freeze-dried powder

**Buffer formulation:** PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.

**Original Concentration:** 500µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 10.7

**Predicted Molecular Mass:** 13.8kDa

**Accurate Molecular Mass:** 18kDa as determined by SDS-PAGE reducing conditions.

### [ USAGE ]

Reconstitute in ddH<sub>2</sub>O to a concentration of 0.1-0.5 mg/mL. Do not vortex.


### [ STORAGE AND STABILITY ]

**Storage:** Avoid repeated freeze/thaw cycles.


Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was

[illegible]

### Figure 1. Gene Sequencing (Extract)



### Figure 2. SDS-PAGE

**[ IMPORTANT NOTE ]**

The kit is designed for in vitro and research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.