

Recombinant Epidermal Growth Factor (EGF)

Catalog No.: TP09104

50µg

Sequence Information

Species: Pig

Gene ID:397083

Swiss Prot:Q00968

Synonyms:URG; Beta-Urogastrone;
Pro-epidermal growth factor

Residues:Ala45~His275

APFLIFSHGNSIFRIDLEGTNHEQLVADAGISVLMDFHYN EERIYWVDLERQLL
QRVFLNGTRQEKVCNLEKNVSGMAINWINEELIWSNQEGTITVTDMKGNNSRV
LLSALKYPANVAVDPVERLMFWSSVVGSLHRADVTGVEVRL LLETSEEIAAVS
LDVLDKRLFWIQYNREGGSSRICSCDYDGGSVHFSKHLTQHNVFAMSLFGDHIF
YSTWKKKTIWVANKH

Product Information

Source: Recombinant expression.

Host: *E.coli*

Tags: N-terminal His-Tag

Subcellular Location: Membrane.

Purity: >90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.8

Predicted Molecular Mass: 33.1kDa

Accurate Molecular Mass: 33kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH₂O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

The kit is designed for in vitro and research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.