

Recombinant Rat Ciliary Neurotrophic Factor (CNTF)

Catalog No.: TP01624

100µg

Sequence Information

Species: Rat

Gene ID:25707

Swiss Prot:P20294

Synonyms:HCNTF

Residues:Ala2~Met200

AFAEQTPLTLHRRDLCRSIWLARKIRSDLTALMESYVKHQGLNKNINLDSVDG

VPVASTDRWSEMTAERLQENLQAYRTFQGMLTKLLEDQRVHFTPTEGDFHQAI

HTLMLQVSAFAYQLEELMVLLEQKIPENEADGMPATVGDGGLFEKKLWGLKVLQ

ELSQWTVRSIHDLRVISSHQMGISALESHYGAKDKQM

Product Information

Source: Recombinant expression.

Host: *E.coli*

Tags: N-terminal His Tag

Subcellular Location: Cytoplasm

Purity: >90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 6.4

Predicted Molecular Mass: 26.2kDa

Accurate Molecular Mass: 28kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH₂O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was

The kit is designed for in vitro and research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.