



Recombinant Human High Temperature Requirement Factor A1 (HTRA1)

Catalog No.: TP09414

100µg

Sequence Information

Species: Human

Gene ID:5654

Swiss Prot:Q92743

Synonyms:L56; ORF480

Residues:Gly204-Leu364

GSGFIVSEDGLIVTNAHVVTNKHVRVKVELKNGATYEAKIKDVDEKADIALIKIDHQGKLPVLLLLGRSSELRPGEFVVAIGSPFSLQNTVTTGIVSTTQRG
GKELGLRNSDMDYIQTDAIINYGNSGGPLVNLDGEVIGINTLKVTAGISFAIPSDKIKKFL

Product Information

Source: Prokaryotic expression.

Host: *E. coli*

Tags: N-terminal His Tags+GST

Subcellular Location: Secreted.

Purity: >90%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% skl,5%Trehalose.

Original Concentration: 50µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 6.5

Predicted Molecular Mass: 47.2kDa

Accurate Molecular Mass: 47kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH₂O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[IDENTIFICATION]




Figure . SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for in vitro and research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.