

Recombinant Odorant Binding Protein 2A (OBP2A)

Catalog No.: TP09824

100µg

Sequence Information

Species: Human

Gene ID:29991

Swiss Prot:Q9NY56

Synonyms:OBP2-A; OBP; OBP2C; OBPIIa;
hOBPIIa

Residues:Leu16~His170

LSFTLEEEDITGTWYVKAMVVDKDFPEDRRPRKVSPVKVTALGGGNLEATFTFM

REDRCIQKKILMRKTEEPGKFSAYGGRKLIYLQELPGTDDYVFYCKDQRRGGLR

YMGKLVGRNPNTNLEALEEFKKLVQHKGLSEEDIFMPLQTGSCVLEH

Product Information

Source: Recombinant expression.

Host: *E.coli*

Tags: N-terminal His and GST Tag

Subcellular Location: Secreted.

Purity: >90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 7.9

Predicted Molecular Mass: 47.8kDa

Accurate Molecular Mass: 50kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH₂O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was

The kit is designed for in vitro and research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.